TDX™

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SPECIFIC INSTRUMENTS

55-1001	Bone awl
55-1002	Bone probe
55-1003	Curved bone probe
55-1004	Straight Sounder
55-1005	Curved Sounder
57-0011	3.5 mm bone tap
57-0010	4.5 mm bone tap
55-1025	5.5 mm bone tap
55-1026	6.5 mm bone tap
55-1027	7.5 mm bone tap
55-1028	8.5 mm bone tap
55-1037	Multi-axial screw driver, wing
55-1060	Set screw driver/holder
55-1062	Set screw Intermediate driver
55-1065	Deflection beam torque wrench (not ratchet)
55-1068	Torque T-Handle 100 in-lbs

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Description: The Orthofix Spinal Fixation System is a temporary, titanium alloy, multiple component system comprised of a variety of non-sterile, single use components that allow the surgeon to build a spinal implant construct. The system is attached to the vertebral body by means of screws, and hooks to the non-cervical spine. The Orthofix Spinal Fixation System consists of an assortment of screws, hooks, rods, spacers, staples, washers, dominos, lateral offsets, and cross connectors. The Orthofix TDX Stabilization System, which can be used with non-reduction style pedicle screws and cross connectors, is a sterile, single-use component that allows the surgeon to build a spinal implant construct. The TDX Stabilization System consists of polyurethane core which is secured in an implant grade titanium housing. The Orthofix Spinal Fixation System titanium implants are not compatible with components or metal from any other manufacturer's system.

Levels of Use: The Orthofix Spinal Fixation System is intended for non-cervical use in the spine. When used as a non-pedicle anterolateral fixation system it may be used from levels T1 to S1. When used with pedicle screw fixation, the Orthofix Spinal Fixation System will be used at L5-S1, with the device fixed or attached to the lumbar and sacral spine (levels of pedicle screw fixation are L3 and below). When used as a posterior non-pedicle fixation system it may be used from levels T1-S1. When used as a non-pedicle anterolateral screw fixation system to the non-cervical spine, the staple and washer may be used from levels T6 to L5.

Indications	Applicability
<p>The Orthofix Spinal Fixation System is intended for non-cervical use in the spine. The Orthofix Spinal Fixation System, when used for pedicle screw fixation, is intended only for patients:</p> <ul style="list-style-type: none"> a) Having severe spondylolisthesis (Grades 3 and 4) at the L5-S1 joint; b) Who are receiving fusion using autogenous bone graft only; c) Who are having the device fixed or attached to the lumbar and sacral spine (L3 and below); and d) Who are having the device removed after the development of a solid fusion mass. 	<p>Pedicle Screws (non-reduction style), Rods, Cross Connectors, Dominos, Lateral Offsets, and Spacers.</p>
<p>The Orthofix Spinal Fixation System, when used as a pedicle screw system in skeletally mature patients, is intended to provide immobilization and stabilization of spinal segments, as an adjunct to fusion in the treatment of the following acute and chronic instabilities or deformities of the thoracic, lumbar and sacral spine:</p> <ul style="list-style-type: none"> a) Degenerative spondylolisthesis with objective evidence of neurologic impairment; b) Fracture; c) Dislocation; d) Scoliosis; e) Kyphosis; f) Spinal tumor; and g) Failed previous fusion (pseudoarthrosis). 	<p>Pedicle Screws (non-reduction style), Rods, Cross Connectors, Dominos, Lateral Offsets, and Spacers.</p>

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<p>The Orthofix Spinal Fixation System, when used for anterolateral non-pedicle screw fixation to the non-cervical spine, is intended for the following indications:</p> <ul style="list-style-type: none"> a) Degenerative disc disease (as defined as back pain of discogenic origin with degenerative disc confirmed by history and radiographic studies); b) Spondylolisthesis; c) Spinal stenosis; d) Spinal deformities (i.e., scoliosis, kyphosis, and/or lordosis); e) Tumor; f) Pseudoarthrosis; g) Failed previous fusion; and h) Trauma (i.e., fracture or dislocation). 	<p>Screws (non-reduction style), Rods, Cross Connectors, Dominos, Lateral Offsets, Spacers, Staples, and Washers.</p>
<p>The Orthofix Spinal Fixation System, when used for posterior non-pedicle screw fixation system of the non-cervical spine, is intended for the following indications:</p> <ul style="list-style-type: none"> a) Degenerative disc disease (as defined as back pain of discogenic origin with degenerative disc confirmed by history and radiographic studies); b) Spondylolisthesis; c) Spinal stenosis; d) Spinal deformities (i.e., scoliosis, kyphosis, lordosis); e) Tumor; f) Pseudoarthrosis; g) Failed previous fusion; and h) Trauma (i.e., fracture or dislocation). 	<p>Hooks, Rods, Cross Connectors, and Dominos.</p>
<p>When used as a pedicle screw fixation system in skeletally mature patients, the TDX Stabilization System is intended to provide immobilization and stabilization of spinal segments as an adjunct to fusion in the treatment of the following acute and chronic instabilities for deformities of the thoracic, lumbar, and sacral spine: degenerative spondylolisthesis with objective evidence of neurologic impairment, kyphosis, and failed previous fusion (pseudoarthrosis).</p> <p>In addition, when used as a pedicle screw fixation system, the TDX Stabilization System is indicated for use in patients:</p> <ul style="list-style-type: none"> • Who are receiving fusions with autogenous graft only; • Who are having the device fixed or attached to the lumbar or sacral spine; and • Who are having the device removed after the development of a solid fusion mass 	<p>TDX Stabilization System</p>

Note: For all of these indications, bone graft must be used.

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Contraindications include, but are not limited to:

1. Morbid obesity
2. Mental Illness
3. Alcoholism or drug abuse
4. Pregnancy
5. Metal sensitivity/allergies
6. Severe osteopenia
7. Patients unwilling or unable to follow post-operative care instructions
8. If implanting the TDX Stabilization System, known allergy to titanium, polyurethane, or ethylene oxide residuals
9. Any circumstances not listed under the heading indications.

Potential Adverse Events

All of the possible adverse events associated with spinal fusion surgery without instrumentation are possible. With instrumentation, a listing of possible adverse events includes, but is not limited to:

1. Device component fracture
2. Loss of fixation
3. Non-union
4. Fracture of the vertebra
5. Neurological injury
6. Vascular or visceral injury
7. Early or late loosening of any or all of the components
8. Disassembly and/or bending of any or all components
9. Foreign body (allergic) reaction to implants, debris, corrosion products, graft material, including metallosis, staining, tumor formation, and/or auto-immune disease
10. Pressure on the skin from component parts in patients with inadequate tissue coverage over the implant possibly causing skin penetration, irritation, and/or pain
11. Post-operative change in spinal curvature, loss of correction, height, and/or reduction
12. Infection
13. Pain, discomfort, or abnormal sensations due to the presence of the device
14. Hemorrhage
15. Cessation of any potential growth of the operated portion of the spine
16. Death

Note: Potential risks identified with the use of the device system may require additional surgery.

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Warnings and Precautions

- 1) The safety and effectiveness of pedicle screw systems have been established only for spinal conditions with significant mechanical instability or deformity requiring fusion with instrumentation. These conditions are: significant mechanical instability or deformity of the thoracic, lumbar, and sacral spine secondary to degenerative spondylolisthesis with objective evidence of neurological impairment, fracture, dislocation, scoliosis, kyphosis, spinal tumor, and failed previous fusion (pseudoarthrosis). The safety and effectiveness of these devices for any other condition are unknown.
- 2) When used as a pedicle screw implant system, the device system is intended only for grade 3 or 4 spondylolisthesis at the fifth lumbar – first sacral (L5-S1) vertebral joint.
- 3) The screws of this device system are not intended for insertion into the pedicles to facilitate spinal fusion above the L5-S1 joint.
- 4) Benefit of spinal fusions utilizing any pedicle screw fixation system has not been adequately established in patients with stable spines.
- 5) Potential risks identified with the use of this device system, which may require additional surgery, include: device component fracture, loss of fixation, non-union, fracture of the vertebra, neurological injury, and vascular or visceral injury.
- 6) Single use only .
- 7) Non-sterile; the screws, hooks, dominos, lateral offsets, spacers, staples, washers, locking nuts, cross connectors, and instruments are sold non-sterile, and therefore must be sterilized before use.
- 8) The TDX Stabilization System is sold sterile, and therefore should not be resterilized.
- 9) To facilitate fusion, a sufficient quantity of autologous bone or other appropriate material should be used.
- 10) Failure to achieve arthrodesis will result in eventual loosening and failure of the device construct.
- 11) Excessive torque applied to the screws may strip the threads in the bone.
- 12) DO NOT REUSE IMPLANTS. Discard used, damaged, or otherwise suspect implants.
- 13) The implantation of pedicle screw spinal systems should be performed only by experienced spinal surgeons with specific training in the use of this pedicle screw spinal system because this is a technically demanding procedure presenting a risk of serious injury to the patient.
- 14) Based on fatigue testing results, the physician/surgeon should consider the levels of implantation, patient weight, patient activity level, other patient conditions, etc. which may impact on the performance of the system.
- 15) The safety and effectiveness of this device for the indication of spinal stabilization without fusion have not been established.
- 16) Reuse of devices labeled as single-use could result in injury or re-operation due to breakage or infection. Do not sterilize single-use implants that come in contact with body fluids.

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