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Referencias

1. **Magnan B et al.** G.I.O.T., 2006; 32 (suppl.1), S322-S326.
2. **Magnan B et al.** J Bone Joint Surg, 2006; 88 B, 1474-9
3. **Cotton FJ, Wilson LT.** Fractures of the os calcis. Boston Medical Journal, 1908; 159, 559-565.
4. **Howard JL, Buckley R, McCormack R, Pate G, Leighton R, Petrie D, Galpin R.** Complications following management of displaced intra-articular calcaneal fractures: a prospective randomized trial comparing open reduction internal fixation with nonoperative management. Journal of Orthopaedic Trauma, 2003; 4, 241-249.
5. **Rajkumar P, Henderson AA.** Surgical treatment of displaced intra-articular fractures of the os-calcis. Foot and Ankle, 2003; 9, 3-6.
6. **Sanders R, Fortin P, DiPasquale T, Walling A.** Operative treatment in 120 displaced intra-articular calcaneal fractures. Results using a prognostic computed tomography scan classification. Clinical Orthopaedics and Related Researches, 1993; 290, 87-95.
7. **Randle JA, Kreder HJ, Stephen D, Williams J, Jaglal S, Hu R.** Should calcaneal fractures be treated surgically? A meta-analysis. Clinical Orthopaedics and Related Researches, 2000; 377, 217-27.
8. **Richardson ML, Van Vu M, Vincent LM, Sangeorzan BJ, Benirschke SK.** CT measurement of the calcaneal varus angle in the normal and fractured hindfoot. Journal of Computer Assisted Tomography, 1992; 16, 261-264.
9. **Rosenberg ZS, Feldman F, Singson RD.** Intra-articular calcaneal fractures: computed tomographic analysis. Skeletal Radiology, 1987; 16,105-113.
10. **Magnan B, Montanari M, Bragantini A, Bartolozzi P.** A system of prognostic evaluation of CT imaging of heel fractures: the Score Analysis Verona (SAVE). Foot Disease, 1995; 1, 19-25.
11. **Csizy M, Buckley R, Tough S, Leighton R, Smith J, McCormack R, Pate G, Petrie D, Galpin R.** Displaced intra-articular calcaneal fractures: variables predicting late subtalar fusion. Journal of Orthopaedic Trauma, 2003; 2, 106-112.
12. **Sanders R.:** Displaced intra-articular fractures of the calcaneus. Journal of Bone and Joint Surgery American Volume, 2000; 82, 225-250.
13. **Coughlin MJ.** Calcaneal fractures in industrial patients. Foot Ankle International, 2000; 21 (11): 896-905.
14. **Myerson M, Quill GE Jr.** Late complications of fractures of the calcaneus. Journal of Bone and Joint Surgery American Volume, 1993; 3, 331-341.
15. **Paley D, Hall H.** Intra-articular fractures of the calcaneus. A critical analysis of results and prognostic factors. Journal of Bone and Joint Surgery American Volume, 1993; 3, 342-354.
16. **Borelli J Jr, Torzilli P.** Effect of impact load on articular cartilage: development of an intra-articular fracture model. Journal of Orthopaedic Trauma, 1997; 11, 319.

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- Barnard L, Odegardd JK.** Conservative approach in the treatment of fractures of the calcaneus. *Journal of Bone and Joint Surgery, American Volume*, 1955; 37-A, 1231-1236.
- Benirschke SK, Kramer PA.** Wound healing complications in closed and open calcaneal fractures. *Journal of Orthopaedic Trauma*, 2004; 18, 1-6.
- Castelli F.** The minimally invasive osteosynthesis of intra-articular calcaneus fractures: early results. Abs. to 7th EFORT CONGRESS, June 2005. P2-504.
- Conn HR:** The treatment of fractures of the os calcis. *Journal of Bone and Joint Surgery American Volume*, 1935; 17, 392-405.
- Dooley P, Buckley R, Tough S, McCormack B, Pate G, Leighton R, Petrie D, Galpin B.** Bilateral calcaneal fractures: operative versus nonoperative treatment. *Foot and Ankle International*, 2004; 2, 47-52.
- Heffernan G, Khan F, Awan N, Riordain CO, Corrigan J. A** comparison of outcome scores in os calcis fractures. *Irish Journal of Medical Science*, 2000; 169(2), 127-128.
- Kenwright J.:** Fractures of the calcaneum, *Journal of Bone and Joint Surgery British Volume*, 1993; 75(2), 176-177.
- Lim EV, Leung JP.** Complications of intraarticular calcaneal fractures. *Clinical Orthopaedics and Related Researches*, 2001; 391, 7-16.
- Loucks C, Buckley R.** Bohler's angle: correlation with outcome in displaced intra-articular calcaneal fractures. *Journal of Orthopaedic Trauma*, 1999; 8, 554-558.
- Magnan B, Caudana R, Campacci A, Barzoi A, Molinaroli F, Zonta L, Martelli S.** Follow-up clinico e radiografico mediante T.C. delle fratture di calcagno trattate con mini FEA. *Chirurgia del piede*, 1992; 16, 145-150.
- Morrey BF, Wiedeman GP Jr.** Complications and long-term results of ankle arthrodesis following trauma. *Journal of Bone and Joint Surgery American Volume*, 1980; 5, 777-784.
- Nogarin L, Magnan B, Bragantini A, Rebeccato A, Montanari M, Schiavon R.** Trattamento a cielo chiuso degli affossamenti talamici con mini-fissatori esterni. Fratture del calcagno. *Progressi in medicina e chirurgia del piede*. Bologna. Aulo Gaggi editore. 1994.
- Paley D, Fischgrund J.** Open reduction and circular external fixation of intraarticular calcaneal fractures. *Clinical Orthopaedics and Related Researches*, 1993; 290, 125-131.
- Pajenda G., Chatwami S., Loidl G., Osterman R., Vecsei V.** Percutaneous screw fixation of the displaced calcaneal fractures. *Book of Abstract SICOT/SIROT 2002*, p 308.
- Rammelt S, Zwipp H.** Calcaneus fractures: facts, controversies and recent developments. *Injury*, 2004; 5, 443-461.
- Sanders R. Hansen S.T. Jr, Mc Reynolds I.S.** Trauma to the calcaneus and its tendon: fractures of the calcaneus. In: *Jahss M.H.: Disorders of the foot and ankle*. 2nd edition, Philadelphia, W.B. Saunders Co., 1991.
- Talarico LM, Vito GR, Zyryanov SY.** Management of displaced intraarticular displaced fractures by using external fixation, minimally invasive open reduction, and early weight-bearing. *Journal of Foot and Ankle Surgery*, 2004; 43(1), 43-50.
- Thornes BS, Collins AL, Timlin M, Corrigan J.** Outcome of calcaneal fractures treated operatively and non-operatively. The effect of litigation on outcomes. *Irish Journal of Medical Science*, 2002; 3, 155-157.
- Yong Soo Choi, Ki Soo Kim.** Minimally invasive open reduction and cross screws fixation for Sanders tipe II fractures of the calcaneus. *Book of Abstract SICOT/SIROT 2002*, p 2020.
- Yong Soo Choi, Ki Soo Kim, Hyeon Gyu Boem.** Biomechanical stability of cross screws fixation fixation for Sanders type II calcaneal fractures. *Book of Abstract SICOT/SIROT 2002*, p 2021.
- Zwipp H, Tscherne H, Thermann H, Weber T.** Osteosynthesis of displaced intraarticular fractures of the calcaneus. Results in 123 cases. *Clinical Orthopaedics and Related Researches*, 1993; 290, 76-86.

Manufactured by: ORTHOFIX Srl
Via Delle Nazioni 9
37012 Bussolengo (Verona)
Italy

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