
Ortholo

POLYMER CERCLAGE SYSTEM

SISTEMA DE BANDA DE CERCLAJE

TÉCNICA QUIRÚRGICA

MBA[®]

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Indicaciones de trauma

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- Fracture Fixation Using Cerclage, Research Applied to Surgery; (2015) By S. M. PERREN et al; AO Foundation, Davos, Switzerland; Acta Chir. orthop. Traum. cech., 82, No. 6 p. 389-397
- Periprosthetic fractures of the femur: the stability of the implant dictates the type of treatment; (2010) By Francesco Bigg et al; Orthopaedics and Traumatology Department, San Martino Hospital, Italy; J Orthopaed Traumatol 11:1-5
- Humeral shaft cerclage wiring: a safe technique to prevent radial nerve injury; (2017) By Stephan Grechenig et al; Clinic of Trauma Surgery, University of Regensburg, Germany; Injury, Int. J. Care Injured 48S5 S12-S14
- Stainless steel wire versus FiberWire suture cerclage fixation to stabilize the humerus in total shoulder arthroplasty; (2014) By Niklas Renner et al; Department of Orthopaedics, University of Zurich, Switzerland; J Shoulder Elbow Surg. 1-7
- Extended trochanteric osteotomy: comparison of 3 modes of fixation: metallic wires, cables, plate, about a series of 157 cases; (2018); Jean Louis Prudhon et al; Tardy Centre osteo articulaire, Echirolles, France; SICOT-J, 4, 21
- Debris From Cobalt-Chrome Cable May Cause Acetabular Loosening; (1992) By Scott S. Kelley et al; University of North California; Clinical Orthopaedics and Related Research - Number 285
- Complications of a Cable Grip System; (1996) By Craig D. Silvertown et al; Department of Orthopaedic Surgery, Rush-Presbyterian St. Luke's Medical Center, Chicago, Illinois; The Journal of Arthroplasty Vol. 11 No. 4
- Periprosthetic femoral fractures around hip arthroplasty: Current concepts in their management; (2009) By Eleftherios TSIRIDIS et al; Academic Department of Trauma and Orthopaedics, University of Leeds, Leeds - UK; Hip International /Vol. 19 no. 2, 75-86
- Patellar Fractures: An Innovative Surgical Technique With Transosseous Suture to Avoid Implant Removal; (2015) By Oscar Buezo et al; Mutua Montan sa, Barcelona, Spain; sri.sagepub.com at Univ California Santa Barbara, Surgical Innovation 1-5
- Technique Using Isoelastic Tension Band for Treatment of Olecranon Fractures; (2015) By R Ajfer RA Danoff JR et al; Am J Orthop. 44 (12):542-546
- Cerclage, evolution and potential of a Cinderella technology. An overview with reference to periprosthetic fractures; (2011) By Perren SM et al; Acta Chir Orthop Traumatol Cech 78(3): 190-199
- Cemented Total Hip Replacement Cable Debris and Acetabular Construct Durability; (2009) By Aaron J. Altenburg, MD, et al.; University of Iowa, Iowa City, The Journal Of Bone And Joint Surgery
- Locking Compression Plate and Cerclage Band for Type B1 Periprosthetic Femoral Fractures Preliminary Results at Average 30-Month Follow-Up; (2011) By Huaming Xue, MD et al; Department of Orthopedic Surgery, Yangpu District Central Hospital, Shanghai, China. The Journal Of Arthroplasty Vol. 26 No. 3
- Intraoperative Periprosthetic Femur Fracture: A Biomechanical Analysis of Cerclage Fixation; (2015) By Nicholas B. Frisch, MD et al, Henry Ford Health System Department of Orthopaedic Surgery, Detroit, Michigan, The Journal Of Arthroplasty
- Failure by Immune Reaction to Metal Debris in Total Joint Arthroplasty of Hip and Knee; (2017) By Boddur C et al, Lenox Hill Hospital, USA, Austin Journal Of Musculoskeletal Disorders
- Biomechanical performance of different cable and wire cerclage configurations; By Mark Lenz et al, (2013); AO Research Institute Davos, Switzerland; INTERNATIONAL ORTHOPAEDICS (SICOT) 37:125-130
- Biomechanical Comparison of Three Sternotomy Closure Techniques: Stati Lateral Distraction at the web site of KINAMED; <http://www.kinamed.com/wordpress/wp-content/uploads/SuperCable-Sternal-Closure-Test-Reports-2016.pdf>
- Biomechanical tests of Ortholox Polymer Cable Knot strength of cerclage bands and wires; By James W. Wilson, (1988); Department of Surgical Sciences, University of Wisconsin; Acta Orthopaedica Scandinavica 59:5, 545-547

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